


Crutches



Unit: mm

	Crutches (2 per Set)		
Product Number	CMS-80L	CMS-80M	CMS-80S
Specification	Size L	Size M	Size S
Length (Total Length)	1,175~1,475 (13-Level)	1,075~1,375 (13-Level)	965~1,115 (7-Level)
Length (Underarm ~ Handgrip)	355~530 (8-Level)	305~455 (7-Level)	255~405 (7-Level)
Weight (Each)	1,000g	960g	940g
Color	Silver/Black		
Diameter (Rubber Tip)	19mm (Ø)		
Manufacturer's Suggested Retail Price	¥11,880 (Non-Taxable)		

© Specifications are subject to change without prior notice for continual product improvement.

Quad Cane



Unit: mm

	Quad Cane		
Product Number	CMS-40L	CMS-40M	CMS-40S
Specification	Size L	Size M	Size S
Length	710~935 (10-Level)		
Base Size	315 (290) × 220	265 (250) × 190	220 (205) × 165
Weight	1,060g	1,020g	980g
Color	Silver/Black		
Diameter (Rubber Tip)	16mm (Ø)		
Manufacturer's Suggested Retail Price	¥8,450 (Non-Taxable)	¥8,190 (Non-Taxable)	¥7,920 (Non-Taxable)

© Specifications are subject to change without prior notice for continual product improvement.

Lofstrand Crutches



Unit: mm

	Lofstrand Crutches		Seven Crutches		
Product Number	CMS-60L	CMS-60S	CMS-70L	CMS-70M	CMS-70S
Specification	Size L	Size S	Size L	Size M	Size S
Length (Grip Height)	645~970	520~770	710~960	535~735	420~545
Length (Cuff Length)	200~300	180~260	200~300	180~260	160~220
Cuff Diameter	90mm (L)	70mm (M)	90mm (L)	70mm (M)	65mm (S)
Weight	630g	510g	595g	530g	480g
Color	Silver/Black		Silver/Black		
Diameter (Rubber Tip)	19 (Ø)		19 (Ø)		
Manufacturer's Suggested Retail Price	¥10,300 (Non-Taxable)		¥11,620 (Non-Taxable)		

© Specifications are subject to change without prior notice for continual product improvement.